



Effective Wound Care

We offer two types of dressing that are used in conjunction with our extriCARE® 2400 and extriCARE® 3600 NPWT pumps.

First our anatomically fitted dressings feature a non-adhesive contact layer which can be applied directly to the wound bed. This layer is centered on an adhesive film which creates a tight seal, and may promote healing. A drainage tube connects each dressing to a canister connected to our pump to remove wound material. They come in a variety of sizes to fit most needs.

Next our standard foam dressings come in three sizes in a kit that includes adhesive film, connector tubing, and a suction bell.

